

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Microtubule-Associated Protein Tau F/MAPT-F/TAU-F Protein

产品货号: TP05528

第三版

别名: Microtubule-associated protein tau; MAPTL; Neurofibrillary tangle protein; MTBT1; Paired helical filament-tau; TAU and MAPT;

描述: Tau proteins are proteins which contain four Tau/MAP repeats. They promote microtubule assembly and stability, and might be involved in the establishment and maintenance of neuronal polarity. They are abundant in neurons of the central nervous system and are less common elsewhere, but are also expressed at very low levels in CNS astrocytes and oligodendrocytes. The tau proteins are the product of alternative splicing from a single gene that in humans is designated MAPT. When tau proteins are defective, and no longer stabilize microtubules properly, they can result in several neurodegenerative disorders such as Alzheimer's disease, Pick's disease, frontotemporal dementia, cortico-basal degeneration and progressive supranuclear palsy.

配方: PBS

分子量: 45.9 kDa

序列: Met 1-Leu441

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.