

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Grancalcin/GCA Protein

产品货号: TP05550

第三版

别名:Grancalcin; GCA; GCL

**描述:**Grancalcin (GCA) is a cytoplasmic granule membrane protein that contains 4 EF-hand domains. GCA is calcium-binding protein and particularly abundant in human neutrophils. GCA is highly expressed in bone marrow, and it can be detected in neutrophils and macrophages. Calcium-binding protein GCA cooperates with SRI and LCP1, so it may play a role in the adhesion of neutrophils to fibronectin. GCA also may play a role in the formation of focal adhesions.

**配方:**PBS

**分子量:**24.0 kDa

**序列:**Met 1-Ile217

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.