



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Chymotrypsin-Like Protease CTRL-1/CTRL Protein(C-6His)

产品货号: TP05554

第三版

别名:Chymotrypsin-Like Protease CTRL-1; CTRL; CTRL1

描述:Chymotrypsin-Like Protease CTRL-1 is a protease that belongs to the peptidase S1 family. Human CTRL-1 is synthesized as a 264 amino acid (aa) precursor that contains an 18 aa signal sequence, 15 aa activation peptide and a 231 aa mature chain. CTRL-1 Contains 1 peptidase S1 domain. It has many molecular functions, such as hydrolase, protease, and serine protease. CTRL-1 plays a role in digest and hydrolyze proteins in biological process.

配方:PBS

分子量:27.2 kDa

序列:Cys19-Asn264

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.