



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human C-X-C Motif Chemokine 14/CXCL14 Protein

产品货号: TP05557

第三版

别名:C-X-C Motif Chemokine 14; Chemokine BRAK<sub>m</sub> MIP-2G; Small-Inducible Cytokine B14; CXCL14; MIP2G; NJAC; SCYB14

描述:Human Chemokine (C-X-C Motif) Ligand 14 (CXCL14) is constitutively expressed in certain normal tissues but is reduced or absent from many established tumor cell lines and human cancers. CXCL14 is known to be a chemoattractant for monocyte and dendritic cells. CXCL14 inhibits angiogenesis and exhibits antimicrobial activities. Mature human and mouse CXCL14 differ by only 2 amino acid residues.

配方:PBS

分子量:9.4 kDa

序列:Ser35-Glu111

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.