

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human C-C Motif Chemokine 27/CCL27 Protein

产品货号: TP05569

第三版

别名:C-C Motif Chemokine 27; CC Chemokine ILC; Cutaneous T-Cell-Attracting Chemokine; CTACK; Eskine; IL-11 R-Alpha-Locus Chemokine; Skinkine; Small-Inducible Cytokine A27; CCL27; ILC; SCYA27

描述:Human Chemokine (C-C Motif) Ligand 27 (CCL27) is a small cytokine that is a member of the CC chemokine family; it is expressed in numerous tissues, including gonads, thymus, placenta and skin. CCL27 elicits its chemotactic effects by binding to the chemokine receptor CCR10. Predominantly expressed in the skin, CCL27 is associated with T cell-mediated inflammation of the skin. Human and Mouse CCL27 share 84% sequence identity in the mature form.

配方:PBS

分子量:10.1 kDa

序列:Phe25-Gly112

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.