

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human C-C Motif Chemokine 3/CCL3 Protein

产品货号: TP05570

第三版

别名:C-C Motif Chemokine 3; G0/G1 Switch Regulatory Protein 19-1; Macrophage Inflammatory Protein 1-Alpha; MIP-1-Alpha; PAT 464.1; SIS-Beta; Small-Inducible Cytokine A3; Tonsillar Lymphocyte LD78 Alpha Protein; CCL3; GOS19-1; MIP1A; SCYA3

描述:Human Chemokine (C-C Motif) Ligand 3 (CCL3) is a small cytokine belonging to the CC chemokine family. CCL3 is primarily expressed in T cells, B cells, and monocytes after antigen or mitogen stimulation. CCL3 exhibits chemoattractive and adhesive effects on lymphocytes. CCL3 exerts multiple effects on hematopoietic precursor cells and inhibits the proliferation of hematopoietic stem cells in vitro as well as in vivo. CCR1 and CCR5 have been identified as functional receptors for CCL3.

配方:PBS

分子量:7.5 kDa

序列:Ser24-Ala92

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.