

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human B2M / Beta-2-microglobulin Protein (His tag)

产品货号: TP05573

第三版

别名:Beta-2 microglobulin

描述:B2M, also known as β 2-Microglobulin or CDABP0092, is a component of MHC class I molecules found expression in all nucleated cells (excludes red blood cells). The major function of MHC class I molecules is to display fragments of proteins from within the cell to T-cells and cells containing foreign proteins will be attacked. B2M(β 2-Microglobulin) is a low molecular weight protein. It was demonstrated that B2M(β 2-Microglobulin) was localized in the membranes of nucleated cells and was found to be associated with HL-A antigens. B2M(β 2-Microglobulin) is present in free form in various body fluids and as a subunit of histocompatibility antigens on cell surfaces lateral to the β 3 chain. Unlike β 3, β 2 has no transmembrane region. Directly above β 2 lies the β 1 chain, which itself is lateral to the β 2. In the absence of B2M(β 2 microglobulin), very limited amounts of MHC class I (classical and non-classical) molecules can be detected on the surface. In the absence of MHC class I, CD8 T cells, a subset of T cells involved in the development of acquired immunity cannot develop. Low levels of B2M(β 2 microglobulin) can indicate non-progression of HIV.

配方:PBS

分子量:13.5 kDa

序列:Met 1-Met 119

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.