



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TXNRD1 / TRXR1 Protein (aa 161-647, His tag)

产品货号: TP05589

第三版

别名:GRIM-12;TR;TR1;TRXR1;TXNR

描述: Thioredoxin reductase 1 (TXNRD1) which is a selenocysteine-containing protein is overexpressed in many malignancies. TXNRD1 plays a key role in regulating cell growth and transformation, and protects cells against oxidative damage. We investigated the association between TXNRD1 polymorphisms and ATDH susceptibility. Moreover, TXNRD1 is an essential selenium-containing enzyme involved in detoxification of reactive oxygen species (ROS) and redox signaling. And genetic variations in TXNRD1 favor the development of Drug-induced liver injury (DILI), which is the most common adverse drug reaction.

配方: PBS

分子量: 55 kDa

序列: Tyr 161-Cys 647

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.