

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ELK1 Protein (His & GST tag)

产品货号: TP05600

第三版

别名:RP23-275N2.2;Elk-1

描述:E twenty-six (ETS)-like transcription factor 1, also known as Elk1 or Member of ETS oncogene family (ELK1), is a member of the ETS oncogene superfamily which is characterized by a common protein domain that regulates DNA binding to target sequences. Elk1 is expressed in the nuclei of non-neuronal cells and function as a transcription activator. It plays important roles in various contexts, including long-term memory formation, drug addiction, Alzheimer's disease, Down syndrome, breast cancer, and depression.

配方:PBS

分子量:73 kDa

序列:Met 1-Pro 428

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.