
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Smad3 Protein (His & GST tag)

产品货号: TP05613

第三版

别名:AU022421;Madh3

描述:SMAD3 belongs to the SMAD family. Members of this family mediate signal transduction by the TGF-beta/activin/BMP-2/4 cytokine superfamily from receptor Ser/Thr protein kinases at the cell surface to the nucleus. SMAD3 is involved in cell signalling. It modulates signals of activin and TGF β 's. Binding of SMAD3 with SMAD4 enables its transmigration into the nucleus where it forms complexes with other proteins and acts as a transcription factor. SMAD3 is a receptor-regulated SMAD (R-SMAD). In mice, mutation of SMAD3 has been linked to colorectal adenocarcinoma, increased systemic inflammation, and accelerated wound healing. Increased SMAD3 activity has been implicated in the pathogenesis of scleroderma. Smad3 is also a multifaceted regulator in adipose physiology and the pathogenesis of obesity and type 2 diabetes.

配方:PBS

分子量:75.9 kDa

序列:Met1-Ser425

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.