

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Protein A/G (Phycoerythrin)

产品货号: TP05665

第三版

别名:

**描述:**The recombinant Protein A/G consists of 7 IgG-binding domains EDABC-C1C3, which corresponds to the Protein A and G domains that are included in the recombinant sequence. The Protein A portion is from Staphylococcus aureus segments E, D, A, B and C. The Protein G portion is from Streptococcus segments C1 and C3. The binding capacity of recombinant Protein A/G is broader than either Protein A or Protein G alone.

**配方:**PBS

**分子量:**

**序列:**

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.