



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Growth Arrest-Specific Protein 7/GAS-7/KIAA0394 Protein(N-6His)

产品货号: TP05669

第三版

别名:Growth Arrest-Specific Protein 7; GAS-7; GAS7; KIAA0394

**描述:**Growth Arrest-Specific Protein 7 (GAS7) is expressed primarily in terminally differentiated brain cells and predominantly in mature cerebellar Purkinje neurons. GAS7 may play a role in neuronal development by promoting maturation and morphological differentiation of cerebellar neurons. Inhibition of GAS7 production in terminally differentiating cultures of embryonic murine cerebellum impedes neurite outgrowth. The hyper-expression of GAS7 may play an major role in the initiation and development of human osteosarcoma.

**配方:**PBS

**分子量:**49.4 kDa

**序列:**Met 1-Ile412

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.