

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse C-X-C motif Chemokine 15/CXCL15/Lungkine Protein(C-6His)

产品货号: TP05680

第三版

别名:C-X-C motif chemokine 15;Cxcl15;Lungkine;Small-inducible cytokine B15;Scyb15

描述:Mouse C-X-C motif chemokine 15, also known as Cxcl15, is a secreted protein which is member of the ELR motif-containing CXC chemokines. It expressed at low levels in fetal lung, the expression restricted to the lung, produced by bronchoepithelial cells and is released into the airways. It plays an important role in lung-specific neutrophil trafficking during normal and inflammatory conditions. It also appears chemotactic for neutrophils.

配方:PBS

分子量:17.4 kDa

序列:Gln26-Ala167

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.