



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human HDGFRP3/HRP-3/HDGF2 Protein(C-6His)

产品货号: TP05685

第三版

别名:Hepatoma-Derived Growth Factor-Related Protein 3; HRP-3; Hepatoma-Derived Growth Factor 2; HDGF-2; HDGFRP3; HDGF2

描述:Hepatoma-Derived Growth Factor-Related Protein 3 (HDGFRP3) belongs to the HDGF family. HDGFRP3 can be found in testis, heart, spinal cord and brain. HDGFRP3 localizes to the nucleus and contains one PWWP domain. HDGFRP3 enhances DNA synthesis and may have a role in cell proliferation.

配方:PBS

分子量:23.6 kDa

序列:Met 1-Thr 203

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.