

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Heat Shock Factor Protein 2/HSF2 Protein(N-6His)

产品货号: TP05690

第三版

别名:Heat Shock Factor Protein 2; HSF 2; Heat Shock Transcription Factor 2; HSTF 2; HSF2; HSTF2

描述:Heat Shock Factor Protein 2 (HSF2) belongs to the HSF family of transcription factors that bind specifically to the heat-shock promoter element and activate transcription. In higher eukaryotes, HSF is unable to bind to the HSE unless the cells are heat shocked. HSF2 is widely expressed in many cells and tissues. HSF2 is located on Cytoplasmic during normal growth. But when it is activated, HSF2 moves to the nucleus.

配方:PBS

分子量:15.9 kDa

序列:Ser411-Ser536

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.