



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Cystatin 8/CST8 Protein(N-6His)

产品货号: TP05691

第三版

别名:Cystatin-8; Cystatin-related epididymal spermatogenic protein; Cystatin-related epididymal-specific protein

描述:Cystatin-8 is a secreted protein which belongs to the cystatin family. It is Lower expression in the testis. Within the testis it is localized to the elongating spermatids, whereas within the epididymis it is exclusively synthesized by the proximal caput epithelium. It performs a specialized role during sperm development and maturation.

配方:PBS

分子量:13.3 kDa

序列:Met1-Phe98

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.