

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Sulfotransferase 4A1/SULT4A1 Protein

产品货号: TP05693

第三版

**别名:**Sulfotransferase 4A1; ST4A1; Brain Sulfotransferase-Like Protein; hBR-STL; hBR-STL-1; Nervous System Sulfotransferase; NST; SULT4A1; SULTX3

**描述:**Sulfotransferase 4A1 (ST4A1) is a member of the Sulfotransferase 1 family. ST4A1 is highly expressed in the cerebral cortex and frontal lobe, but no expression is detected in the pancreas. ST4A1 is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. ST4A1 acts on catecholamines and T4 in a manner that may not involve sulfonation. ST4A1 may have a role in the metabolism of drugs and neurotransmitters in the CNS. In addition, ST4A1 is related to schizophrenia.

**配方:**PBS

**分子量:**33.0 kDa

**序列:**Met 1-Leu284

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.