



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Neurotrimin/NTRI Protein(C-6His)

产品货号: TP05695

第三版

别名: Neurotrimin; NTM; IgLON family member 2; IGLON2; HNT; NTRI;

描述: Neurotrimin is a glycoprotein. It belongs to the immunoglobulin superfamily and IgLON family. Neurotrimin localizes to the cell membrane and contains three Ig-like C2-type domains. Neurotrimin acts as a glycosylphosphatidylinositol (GPI)-anchored neural cell adhesion molecule. Neurotrimin may promote neurite outgrowth and adhesion via homophilic and heterophilic interactions.

配方: PBS

分子量: 32.3 kDa

序列: Gly34-Leu316

纯度: > 95% by HPLC

浓度:

内毒素: < 1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.