

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Hematopoietic Lineage Cell-Specific Protein/HCLS1 Protein(C-6His)

产品货号: TP05704

第三版

别名:Hematopoietic lineage cell-specific protein;Hematopoietic cell-specific LYN substrate 1;LckBP1;p75

描述:Hematopoietic lineage cell-specific protein (HCLS1) is a protein that in humans is encoded by the HCLS1 gene. It is expressed only in tissues and cells of hematopoietic origin. It is substrate of the antigen receptor-coupled tyrosine kinase and plays a role in antigen receptor signaling for both clonal expansion and deletion in lymphoid cells. It may also be involved in the regulation of gene expression.

配方:PBS

分子量:55.0 kDa

序列:Met 1-Glu486

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.