

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Syntaxin-6/STX6 Protein

产品货号: TP05705

第三版

别名:Syntaxin-6; STX6

描述:Syntaxin-6 (STX6) is a single-pass type IV membrane protein, which belongs to the Syntaxin family. STX6 is mainly localized in the plasma membrane. STX6 contains one t-SNARE coiled-coil homology domain and involved in intracellular vesicle trafficking. When STX6 function is inhibited, internalization through caveolae is dramatically reduced, whereas other endocytic mechanisms are unaffected. It is reported that STX6 is necessary for proper expression of focal adhesion kinase and integrin alpha5 adhesion receptor.

配方:PBS

分子量:26.8 kDa

序列:Ser2-Gln234

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.