

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human/Mouse/Rat BMP-2 Protein

产品货号: TP05710

第三版

别名: Bone Morphogenetic Protein 2; BMP-2; Bone Morphogenetic Protein 2A; BMP-2A; BMP2; BMP2A

描述: Bone Morphogenetic Protein-2 (BMP-2) is one of the bone-growth regulatory factors that belong to the transforming growth factor-beta (TGF-beta) superfamily of proteins. BMPs are synthesized as large precursor molecules, which are cleaved by proteolytic enzymes. The active form of BMP-2 can consist of a dimer of two identical proteins or a heterodimer of two related bone morphogenetic proteins.

配方: PBS

分子量: 13.3 kDa

序列: Gln283-Arg396

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.