



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TGOLN2/TGN38 Homolog/TGN46 Protein(C-6His)

产品货号: TP05720

第三版

别名:Trans-Golgi network integral membrane protein 2;TGN38 homolog; TGN46; TGN48; Trans-Golgi network protein TGN51;TGOLN2;TGN46; TGN51

描述:This protein may be involved in regulating membrane traffic to and from trans-Golgi network. Isoform TGN46 is widely expressed. Isoform TGN51 is more abundant in fetal lung and kidney. Isoform TGN48 is barely expressed in embryonic kidney and promyelocytic cells.

配方:PBS

分子量:38.3 kDa

序列:Ala22-Glu381

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.