

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Thimet Oligopeptidase/THOP1 Protein(C-6His)

产品货号: TP05722

第三版

别名:Thimet Oligopeptidase; Endopeptidase 24.15; MP78; THOP1

描述:Thimet Oligopeptidase (THOP1) belongs to the peptidase M3 family which includes neurolysin and mitochondrial intermediate peptidase. THOP1 is located in Cytoplasm. THOP1 is widely expressed in human tissues and can be detected in different subcellular locations. THOP1 is preferential cleavage for bonds with hydrophobic residues at P1, P2 and P3' and a small residue at P1' in substrates of 5 to 15 residues. THOP1 is involved in the metabolism of neuropeptides under 20 amino acid residues and degradation of cytoplasmic peptide. In addition, THOP1 also can degrade the beta-amyloid precursor protein and generate amyloidogenic fragments.

配方:PBS

分子量:80.0 kDa

序列:Lys2-Cys689

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.