



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Clusterin-Like Protein 1/CLUL1 Protein(C-6His)

产品货号: TP05729

第三版

别名:Clusterin-Like Protein 1; Retinal-Specific Clusterin-Like Protein; CLUL1

描述:Clusterin-Like Protein 1 (CLUL1) is a secreted protein that belongs to the clusterin family. CLUL1 is synthesized as a 466 amino acid precursor that contains a 20 amino acid signal sequence, and a 446 amino acid mature chain. CLUL1 is expressed predominantly by cone photoreceptors of the retina. It has been shown that CLUL1 expression is down-regulated in some forms of retinal disease.

配方:PBS

分子量:50.1 kDa

序列:Ala21-Ala442

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.