

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human OX40 Ligand/TNFSF4/OX40L Protein(N-6His)

产品货号: TP05762

第三版

**别名:**Tumor necrosis factor ligand superfamily member 4;Glycoprotein Gp34;OX40 ligand;OX40L;TAX transcriptionally-activated glycoprotein 1;TNFSF4;CD252;TXGP1

**描述:**Tumor necrosis factor ligand superfamily member 4(TNFSF4/OX40L) is a single-pass type II membrane protein.OX40L is expressed on the surface of activated B cells, T cells, dendritic cells and endothelial cells. OX40L binds to OX40 (CD134), a member of the TNF receptor superfamily that is expressed predominantly on activated CD4+ T cells. OX40-OX40L co-stimulates signal to promote the survival and proliferation of activated CD4+ T cells and prolong the immune response. It involved in T-cell proliferation and cytokine production. Additional, it has been found association with systemic lupus erythematosus, no association with occurrence of atherosclerosis.

**配方:**PBS

**分子量:**16.3 kDa

**序列:**Gln51-Leu183

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.