



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PACSIN1/Syndapin-1 Protein(C-6His)

产品货号: TP05764

第三版

别名:Protein Kinase C and Casein Kinase Substrate in Neurons Protein 1; PACSIN1; KIAA1379

描述:Protein Kinase C and Casein Kinase Substrate in Neurons Protein 1 (PACSIN1) belongs to the PACSIN family. PACSIN1 contains one FCH domain and one SH3 domain. PACSIN1 is highly expressed in the brain and at lower levels in the heart, pancreas, and liver. PACSIN1 may play a role in vesicle formation and transport. PACSIN1 has been shown to interact with DNM1, PACSIN3, Huntingtin, and PACSIN2. In addition, PACSIN1 is phosphorylated by casein kinase 2 (CK2) and protein kinase C (PKC).

配方:PBS

分子量:52.0 kDa

序列:Met 1-Ile444

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.