



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PACSIN2/Syndapin-2 Protein(C-6His)

产品货号: TP05765

第三版

别名:Protein Kinase C and Casein Kinase Substrate in Neurons Protein 2; PACSIN2

描述:Protein Kinase C and Casein Kinase Substrate in Neurons Protein 2 (PACSIN2) is a member of the PACSIN family. PACSIN2 is localized to the plasma membrane via its coiled-coil domain. PACSIN2 is widely expressed and contains one FCH domain and one SH3 domain. PACSIN2 forms homo- and hetero-aggregates with other PACSINs. PACSIN2 may play a role in vesicle formation and transport. In addition, PACSIN2 is involved in linking the actin cytoskeleton with vesicle formation by regulating tubulin polymerization.

配方:PBS

分子量:56.7 kDa

序列:Met 1-Gln486

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.