



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Pancreatic Lipase-Related Protein 1/PLRP1 Protein(C-6His)

产品货号: TP05771

第三版

别名:Pancreatic Lipase-Related Protein 1; PL-RP1; PNLIPRP1; PLRP1

描述:Pancreatic Lipase-Related Protein 1 (PNLIPRP1) belongs to the Lipase family within the AB hydrolase superfamily. PNLIPRP1 is a secreted protein and contains one PLAT domain. PNLIPRP1 is involved in lipid metabolic process, acting as a negative regulator of pancreatic lipase activity by competing with pancreatic lipase for colipase occupancy. PNLIPRP1 may play a role in inhibiting dietary triglyceride digestion, but it lacks detectable lipase activity towards triglycerides, diglycerides, phosphatidylcholine, galactolipids or cholesterol esters.

配方:PBS

分子量:51.0 kDa

序列:Lys18-Cys467

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.