

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Transcription Initiation Factor IIB/TFIIB/GTF2B Protein(N-GST)

产品货号: TP05775

第三版

别名:Transcription Initiation Factor IIB; General Transcription Factor TFIIB; S300-II; GTF2B; TF2B; TFIIB

描述:Transcription Initiation Factor IIB (TFIIB) is an essential factor for transcription by RNA Polymerase II. TFIIB localizes to the nucleus where it forms a complex (the DAB complex) with transcription factor IID and IIA. TFIIB plays a role as a bridge between IID, which initially recognizes the promoter sequence, and RNA polymerase II. TFIIB is involved in the selection of transcription start site.

配方:PBS

分子量:61.6 kDa

序列:Met 1-Leu316

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.