

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Parathyroid Hormone-Related Protein/PTHrP Protein(C-6His)

产品货号: TP05784

第三版

别名: Parathyroid Hormone-Related Protein; PTH-rP; PTHrP; Parathyroid Hormone-Like Protein; PLP; PTHrP [1-36]; PTHrP [38-94]; Osteostatin; PTHrP [107-139]; PTHLH; PTHRP

描述: Parathyroid Hormone-Related protein (PTHrP) is a member of the parathyroid hormone family. PTHrP is known as a potent inhibitor of chondrocyte maturation. PTHrP is a critical regulator of cellular and organ growth, development, migration, differentiation and survival and of epithelial calcium ion transport. PTHrP also regulates epithelial-mesenchymal interactions during the formation of the mammary glands. During endochondral bone development, PTHrP plays a major role in suppressing the onset of hypertrophic differentiation. Defects in PTHrP are the cause of Brachydactyly Type E2.

配方: PBS

分子量:

序列: Ala37-Arg175

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.