

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TREM-1/CD354 Protein(C-6His)

产品货号: TP05791

第三版

别名: Triggering Receptor Expressed on Myeloid Cells 1; TREM-1; Triggering Receptor Expressed on Monocytes 1; CD354; TREM1

描述: Triggering Receptor Expressed on Myeloid Cells 1 (TREM-1) is a transmembrane protein with a single Ig-like domain. TREM-1 associates with the adapter protein, DAP12, to deliver an activating signal. TREM-1 is expressed on blood neutrophils and monocytes, and the expression is up-regulated by bacterial LPS. TREM-1 is expressed at high levels on neutrophils of patients with microbial sepsis and in mice with a TREM-1/Fc fusion protein protected mice against LPS-induced shock. Human TREM-1 shares 42% sequence homology with mouse TREM-1.

配方: PBS

分子量: 21.3 kDa

序列: Ala21-Arg200

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.