

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PBLD/MAWBP Protein(N-6His)

产品货号: TP05793

第三版

别名:Phenazine Biosynthesis-Like Domain-Containing Protein; MAWD-Binding Protein; MAWDBP; Unknown Protein 32 From 2D-PAGE of Liver Tissue; PBLD; MAWBP

描述:Phenazine Biosynthesis-Like Domain-Containing protein (PBLD) belongs to the phenazine biosynthesis-like protein (PhzF) family, which is expressed in most tissues. PBLD takes part in the MAPK signaling pathway, and is involved in multiple basic cellular functions. The expression of PBLD can be increased in several disease processes, including insulin resistance, folate deficiency and hypotension.

配方:PBS

分子量:34.0 kDa

序列:Met 1-Ala288

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.