
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse IL-9 Protein

产品货号: TP05795

第三版

别名: Interleukin-9; IL-9; Cytokine P40; T-Cell Growth Factor P40; IL9

描述: Interleukin-9 (IL-9) is a secreted protein that belongs to the IL-7/IL-9 family. Mature mouse IL-9 shares 57% and 74% amino acid sequence identity with human and rat IL-9, respectively. IL-9 supports IL-2 independent and IL-4 independent growth of helper T-cells. IL-9 stimulates cell proliferation and prevents apoptosis. It functions through the IL-9 receptor (IL-9R), which activates different signal transducer and activator (STAT) proteins and thus connects this cytokine to various biological processes. IL-9 has been identified as a candidate gene for asthma. IL-9 is a determining factor in the pathogenesis of bronchial hyperresponsiveness.

配方: PBS

分子量: 15.2 kDa

序列: Gln19-Pro144

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.