



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PDCD10/TFAR15 Protein

产品货号: TP05798

第三版

别名: Programmed Cell Death Protein 10; Cerebral Cavernous Malformations 3 Protein; TF-1 Cell Apoptosis-Related Protein 15; PDCD10; CCM3; TFAR15

描述: Programmed Cell Death Protein 10 (PDCD10) belongs to the PDCD10 family. PDCD10 exists as a homodimer and is widely expressed. PDCD10 can increase mitogen-activated protein kinase activity and MST4 activity. PDCD10 is required for normal cardiovascular development and normal angiogenesis, vasculogenesis and hematopoiesis during embryonic development. Defects in PDCD10 are the cause of cerebral cavernous malformations type 3.

配方: PBS

分子量: 24.9 kDa

序列: Met 1-Ala212

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.