

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human Tryptase alpha/beta-1/TPSAB1 Protein(C-6His)

产品货号: TP05807

---

第三版

别名:trypase alpha/beta 1;Tryptase 1;TPS1; TPS2; TPSAB1; TPSB1

**描述:**Tryptases are serine proteases with trypsin-like specificity. Together with chymases and Cathepsin G, tryptases are important players in mast cell mediation of inflammatory and allergic responses. Tryptase alpha/beta-1(TPSAB1), also known as mast cell protease 7 (MCPT7), it exhibits anticoagulant activity due to its ability to degrade fibrinogen in the presence of a diverse array of protease inhibitors in plasma. The two Isoform 1 and isoform 2 are expressed in lung, stomach, spleen, heart and skin; in these tissues, isoform 1 is predominant. Isoform 2 is expressed in aorta, spleen, and breast tumor, with highest levels in the endothelial cells of some blood vessels surrounding the aorta, as well as those surrounding the tumor and low levels, if any, in mast cells. Isoform 2 cleaves large substrates, such as fibronectin, more efficiently than isoform 1, but seems less efficient toward small substrates. It may play a role in innate immunity.

**配方:**PBS

**分子量:**28.8 kDa

**序列:**Ile31-Pro275

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.