

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human 3-Hydroxybutyrate Dehydrogenase Type 2/BDH2/DHRS6 Protein(N-6His)

产品货号: TP05830

第三版

**别名:**3-Hydroxybutyrate Dehydrogenase Type 2; Dehydrogenase/Reductase SDR Family Member 6; Oxidoreductase UCPA; R-Beta-Hydroxybutyrate Dehydrogenase; BDH2; DHRS6

**描述:**3-Hydroxybutyrate Dehydrogenase Type 2 belongs to the short-chain dehydrogenases/reductases (SDR) family. 3-Hydroxybutyrate Dehydrogenase Type 2 may play an important role in the peripheral utilization of 3-hydroxybutyrate. The cytoplasmic localization with its high ratio of oxidized NAD<sup>+</sup>, the NAD<sup>+</sup> dependence and the kinetic parameters of 3-Hydroxybutyrate Dehydrogenase Type 2 make it suitable to convert high levels of circulating 3-hydroxybutyrate into acetoacetate.

**配方:**PBS

**分子量:**28.9 kDa

**序列:**Met 1-Leu245

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.