

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human 40S Ribosomal Protein S7/RPS7 Protein(C-6His)

产品货号: TP05834

第三版

别名:40S ribosomal protein S7;RPS7

描述:40S ribosomal protein S7(RPS7) belongs to the S7E family of ribosomal proteins. It is phosphorylated by NEK6 during post-translational modification. RPS7 is located in the cytoplasm, binds IPO9 with high affinity. It also can interact with NEK6. As is required for rRNA maturation and typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. The unnormal expression of RPS7 may cause Diamond-Blackfan anemia 8.

配方:PBS

分子量:23.2 kDa

序列:Met 1-Leu194

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.