



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Interleukin-21/IL-21 Protein

产品货号: TP05843

第三版

别名: Interleukin-21; IL21

描述: Interleukin-21 also known as IL-21 is a protein that in mouse is encoded by the IL21 gene, belongs to the IL-15/IL-21 family. Interleukin-21 is a cytokine that has potent regulatory effects on cells of the immune system, including natural killer (NK) cells and cytotoxic T cells that can destroy virally infected or cancerous cells. This cytokine induces cell division/proliferation in its target cells.

配方: PBS

分子量: 14.4 kDa

序列: Pro25-Ser146

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.