

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Peroxiredoxin-6/PRDX6 Protein(N-6His)

产品货号: TP05844

第三版

别名: Peroxiredoxin-6; 1-Cys Peroxiredoxin; 1-Cys PRX; 24 kDa Protein; Acidic Calcium-Independent Phospholipase A2; aiPLA2; Antioxidant Protein 2; Liver 2D PAGE Spot 40; Non-Selenium Glutathione Peroxidase; NSGPx; Red Blood Cells Page Spot 12; PRDX6; AOP2; KIAA0106

描述: Peroxiredoxin 6 belongs to the thiol-specific antioxidant protein family, which is a bifunctional enzyme with two distinct active sites. Peroxiredoxin 6 is involved in redox regulation of the cell and can reduce Hydrogen peroxide and short chain organic, fatty acid, and phospholipid hydroperoxides. Peroxiredoxin 6 may regulates phospholipid turnover and protects against oxidative injury. Peroxiredoxin 6 eases the oxidative stress and TGF-beta-induced abnormalities of human trabecular meshwork cells, it is necessary for blood vessel integrity in injured skin.

配方: PBS

分子量: 27.2 kDa

序列: Met 1-Pro224

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.