

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cystatin S/CST4 Protein(C-6His)

产品货号: TP05847

第三版

别名:Cystatin-S; Cystatin-4; Cystatin-SA-III; Salivary Acidic Protein 1; CST4

描述:Cystatin-S, also called Cystatin-SA-III or Salivary acidic protein, is usually expressed in submandibular and sublingual saliva but not in parotid saliva(at protein level). This protein, which is phosphorylated at both its N- and C-terminal regions, belongs to the cystatin family. Cystatin-S can function as an inhibitor to papain or ficin, it partially inhibits stem bromelain and bovine cathepsin C.

配方:PBS

分子量:15.2 kDa

序列:Ser21-Ala141

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.