



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cytochrome b5 B/CYB5B Protein(C-6His)

产品货号: TP05861

第三版

别名: Cytochrome b5 outer mitochondrial membrane isoform; OMB5

描述: Cytochrome b5 type B (CYB5B) is a membrane of the cytochrome b5 family. It contains 1 cytochrome b5 heme-binding domain. Cytochrome b5 is a membrane bound hemoprotein which function as an electron carrier for several membrane bound oxygenases. In the mitochondrion of eukaryotes and in aerobic prokaryotes, cytochrome b is a component of respiratory chain complex III also known as the bc1 complex or ubiquinol-cytochrome c reductase.

配方: PBS

分子量: 14.0 kDa

序列: Lys12-Cys118

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.