

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ubiquitin-Associated Protein 1/UBAP1/NAG20 Protein(C-6His)

产品货号: TP05867

第三版

别名: Ubiquitin-Associated Protein 1; UBAP-1; Nasopharyngeal Carcinoma-Associated Gene 20 Protein; UBAP1

描述: Ubiquitin-Associated Protein 1 (UBAP1) belongs to the UBA domain family. Members of this family are related to ubiquitin and the ubiquitination pathway. Because of their cytogenetic location, this UBA domain family member is being studied as a putative target for mutation in nasopharyngeal carcinomas. UBAP1 is highly expressed in the heart, brain, placenta, lung, skeletal muscle, liver, and pancreas. UBAP1 consists of two UBA domains and one UMA domain. The ubiquitin associated domain is thought to be a non-covalent ubiquitin binding domain, including a compact three helix bundle.

配方: PBS

分子量: 56.2 kDa

序列: Met 1-Ser502

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.