

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cytochrome b-c1 Complex Subunit 6/UQCRH Protein(N-GST)

产品货号: TP05868

第三版

别名:Cytochrome b-c1 complex subunit 6; mitochondrial;Complex III subunit 6;Complex III subunit VIII;Cytochrome c1 non-heme 11 kDa protein;Mitochondrial hinge protein;Ubiquinol-cytochrome c reductase complex 11 kDa protein

描述:Cytochrome b-c1 complex subunit 6, mitochondrial (UQCRH) belongs to the UQCRH/QCR6 family, it is a subunit of the respiratory chain protein Ubiquinol Cytochrome c Reductase. This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. UQCRH may mediate formation of the complex between cytochromes c and c1. It may mediate formation of the complex between cytochromes c and c1.

配方:PBS

分子量:37.5 kDa

序列:Met 1-Lys91

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.