



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cytoplasmic Dynein Light Chain 1/DYNLL1 Protein(N-6His)

产品货号: TP05875

第三版

别名: Dynein Light Chain 1 Cytoplasmic; 8 kDa Dynein Light Chain; DLC8; Dynein Light Chain LC8-Type 1; Protein Inhibitor of Neuronal Nitric Oxide Synthase; PIN; DYNLL1; DLC1; DNCL1; DNCLC1; HDLC1

描述: Human Dynein Cytoplasmic Light Chain 1 (DYNLL1) has been identified as a protein that interacts with NOS1, leading to NOS1 inhibition. NOS1 dimer is destabilized after binding DYNLL1 a conformation necessary activity, and it regulate numerous biologic processes through its effects on nitric oxide synthase activity. DYNLL1 is widely expressed, with higher expression in testis and moderate expression in brain.

配方: PBS

分子量: 12.5 kDa

序列: Met 1-Gly89

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.