



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Placental Lactogen/CSH1 Protein(C-6His)

产品货号: TP05878

第三版

别名:Chorionic Somatomammotropin Hormone; Choriomammotropin; Lactogen; Placental Lactogen; PL; CSH1; CSH2

描述:Chorionic Somatomammotropin Hormone (CSH1) belongs to the Somatotropin/Prolactin family. It is located at the growth hormone locus on chromosome 17. It is produced by cells in the syncytiotrophoblast layer of the placenta only during pregnancy. Although CSH1 does not interact with GHR but only activates PRLR through zinc-induced dimerization. It is expressed mainly in the placenta and utilizes multiple transcription initiation sites. CSH1 plays an important role in stimulating lactation, growth control and metabolism.

配方:PBS

分子量:23.4 kDa

序列:Val27-Phe217

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.