

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Recombinant Human Ubiquitin-Conjugating Enzyme E2 J2/UBE2J2 Protein(N-GST)

产品货号: TP05879

第三版

别名: Ubiquitin-Conjugating Enzyme E2 J2; Non-Canonical Ubiquitin-Conjugating Enzyme 2; NCUBE-2; UBE2J2; NCUBE2

描述: Ubiquitin-Conjugating Enzyme E2 J2 (UBE2J2) belongs to the ubiquitin-conjugating enzyme family. UBE2J2 is involved in the ubiquitination. UBE2J2 located in the membrane of the endoplasmic reticulum, catalyzes the covalent attachment of ubiquitin to other proteins. UBE2J2 may play a important role in the selective degradation of misfolded membrane protein from the endoplasmic reticulum.

配方: PBS

分子量: 38.0 kDa

序列: Glu124-Arg224

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.