

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human DCN1-Like Protein 1/DCUN1D1 Protein(N-6His)

产品货号: TP05881

第三版

别名:DCN1-Like Protein 1; DCUN1 Domain-Containing Protein 1; Defective in Cullin Neddylation Protein 1-Like Protein 1; Squamous Cell Carcinoma-Related Oncogene; DCUN1D1; DCUN1L1; RP42; SCCRO

描述:DCN1-Like Protein 1 is a protein containing 1 DCUN1 domain and 1 UBA-like domain. DCN1-Like Protein 1 contains part of an E3 Ubiquitin Ligase Complex for Neddylation. It is required for Neddylation of Cullin components of E3 Cullin-RING Ubiquitin Ligase complexes. It enhances the rate of Cullins Neddylation. DCN1-Like Protein 1 recruits the NEDD8-charged E2 Enzyme to the Cullin component. DCUN1D1 is involved in the release of inhibitory effects of CAND1 on the Cullin-RING Ligase E3 complex assembly and activity. It acts also as an oncogene facilitating malignant transformation and carcinogenic progression.

配方:PBS

分子量:32.3 kDa

序列:Met 1-Val259

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.