

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Platelet-Derived Growth Factor BB/PDGF-BB Protein

产品货号: TP05888

第三版

别名:Platelet-Derived Growth Factor Subunit B; PDGF Subunit B; PDGF-2; Platelet-Derived Growth Factor B Chain; Platelet-Derived Growth Factor Beta Polypeptide; Proto-Oncogene c-Sis; Becaplermin; PDGFB; PDGF2; SIS

描述:Platelet-Derived Growth Factor Subunit B (PDGFB) belongs to the PDGF/VEGF growth factor family. Platelet-derived growth factor is a potent mitogen for cells of mesenchymal origin. PDGFB can exist either as a homodimer (PDGF-BB) or as a heterodimer with the platelet-derived growth factor alpha polypeptide (PDGF-AB), where the dimers are connected by disulfide bonds. Mutations in this gene are associated with meningioma. Binding of PDGFB to its receptor elicits a variety of cellular responses. In addition, PDGFB is released by platelets upon wounding and plays an important role in stimulating adjacent cells to grow and thereby heals the wound.

配方:PBS

分子量:12.4 kDa

序列:Ser82-Thr190

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.