

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Interleukin-21/IL-21 Protein

产品货号: TP05892

第三版

别名: Interleukin-21; IL-21; Za11; IL21

描述: IL-21 induces the production of IgG1 and IgG3 in B-cells. IL-21 may promote the transition between innate and adaptive immunity. IL-21 may have a role in proliferation and maturation of natural killer cells in synergy with IL15. IL-21 stimulates interferon gamma production in T-cells and natural killer cells by synergy with IL15 and IL18. Dysregulation of IL-21 has a role in multiple immune-mediated diseases including lupus, psoriasis and chronic inflammatory diseases.

配方: PBS

分子量: 15.6 kDa

序列: Gln23-Ser155

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.