



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Interleukin-28B/IL-28B/IFN-lambda 3 Protein(C-6His)

产品货号: TP05904

第三版

别名: Interleukin-28B; IL-28B; Cytokine Zcyto22; Interferon Lambda-3; IFN-Lambda-3; Interferon Lambda-4; IFN-Lambda-4; Interleukin-28C; IL-28C; IL28B; IFNL3; IFNL4; IL28C; ZCYTO22

描述: Interleukin-28B, also known as Cytokine Zcyto22, Interferon lambda-3, Interferon lambda-4, IFNL3, IFNL4, ZCYTO22 and IL28B, is a secreted cytokine which belongs to the IL-28/IL-29 family. IL-28 has also been shown to play a role in the adaptive immune response. IL28B has immunomodulatory activity and up-regulates MHC class I antigen expression. IL28B displays potent antiviral activity and antitumor activity. In addition, IL28B is a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IL28RA. The ligand/receptor complex seems to signal through the Jak-STAT pathway.

配方: PBS

分子量: 20.7 kDa

序列: Val22-Val196

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.